

# *Fusion QbD<sup>®</sup> – DOE Experimental Design*



S-Matrix Corporation  
1594 Myrtle Avenue  
Eureka, CA 95501  
USA  
Phone: 707-441-0404  
URL: [www.smatrix.com](http://www.smatrix.com)

# Regulatory Guidances – DOE is a Central Element

ICH Q8(R2) Guidance for Industry, Pharmaceutical Development, August, 2009

## Formal Experimental Design:

A structured, organized method for **determining the relationship between factors** affecting a *process* and the *output of that process*. Also known as “Design of Experiments”.

Check out the recent article titled “A Design for Life Sciences” using the link below.  
(The Column, Volume 15, Issue 6, pg 2–7).

<http://www.chromatographyonline.com/design-life-sciences>

# Regulatory Guidances – DOE is a Central Element

Proposed New USP General Chapter: The Analytical Procedure Lifecycle (1220)

## Design of experiments (DOE)

is a **fundamental methodology** for the QRM process. It is a systematic method to determine the relationships between variables affecting a process, and it is used to find cause-and-effect relationships.

...

DOE also utilizes statistical data treatment, which allows clear determinations regarding the significance of a variable and/or its interactions towards the output.

## Supports All These Separation Modes

Reversed Phase

Normal Phase


Chiral

HILIC

Ion Exchange

Size Exclusion

SFC



# Fusion QbD – DOE Study Parameter Flexibility

QbD – requires ability to study any combination of LC parameters which can interactively effect your performance requirements:


- Isocratic Methods
- Gradient Methods
- Any pump program steps – e.g.
  - Equilibration Time & %
  - Isocratic Hold Time & %
  - Gradient Time & Slope
  - Final Hold Time & %
  - Re-equilibration Time & %
- Strong Solvent Type
- Mobile Phase Solvent Blend
- Column Temperature
- pH
- Column Type
- Flow Rate
- Injection Volume
- Ionic Strength

# DOE = Statistical Sampling (Efficient Multifactor Studies)

Consider two variables – five study levels each:


**Brute Force Approach =  
All Possible Combinations**

25 Methods




**Design of Experiments =  
Statistical Sampling**

13 Methods




# DOE – Adding a 3rd and 4th Factor to the Study

3<sup>rd</sup> Variable – pH




**Studies can become very complex – difficult to carry out**


4<sup>th</sup> Variable – Column Type



C18



Phenyl



C8

# Fusion QbD – Experiment Automation to the Rescue!

Built-in Design Wizard which **Automatically Selects the most efficient DOE design** for your stage of work (e.g. chemistry screening or method optimization) and the LC parameters you selected for study.


- **Automatically organizes the design for efficient execution.**
- **Automatically builds the instrument methods and sequence – automated error-free execution.**

	Run No.	Gradient Time	pH	Column Type
1	Condition Column - 1	2	2.73	BEH C18
2	Condition Column - 2	2	2.73	BEH Shield RP18
3	Condition Column - 3	2	2.73	HSS T3
4	Condition Column - 4	2	2.73	CSH Phenyl Hexyl
5	1	25	2.73	BEH C18
6	2	10	2.73	BEH C18
7	3	25	2.73	BEH Shield RP18
8	4	10	2.73	BEH Shield RP18
9	5	25	2.73	HSS T3
10	6	10	2.73	HSS T3
11	7	17.5	2.73	CSH Phenyl Hexyl
12	8	17.5	2.73	CSH Phenyl Hexyl
13	Condition Column - 5	2	3.2	BEH C18
14	Condition Column - 6	2	3.2	BEH Shield RP18
15	Condition Column - 7	2	3.2	HSS T3
16	Condition Column - 8	2	3.2	CSH Phenyl Hexyl
17	9	21.3	3.2	BEH C18
18	10	13.8	3.2	HSS T3
19	11	17.5	3.2	BEH Shield RP18
20	12	17.5	3.2	CSH Phenyl Hexyl
21	Condition Column - 9	2	3.69	BEH C18
22	Condition Column - 10	2	3.69	BEH Shield RP18
23	Condition Column - 11	2	3.69	HSS T3
24	Condition Column - 12	2	3.69	CSH Phenyl Hexyl
25	13	17.5	3.69	BEH C18
26	14	17.5	3.69	BEH Shield RP18
27	15	17.5	3.69	HSS T3
28	16	17.5	3.69	CSH Phenyl Hexyl
29	17	25	3.69	CSH Phenyl Hexyl
30	18	10	3.69	CSH Phenyl Hexyl
31	19	17.5	3.69	BEH C18
32	20	17.5	3.69	BEH Shield RP18

**Auto-built for walk-away execution in the CDS.**



# Fusion QbD – Automated DOE Experiment Workflow




Fusion QbD



Generates QbD-aligned  
DOE Experiment



Automatically Builds  
Sequence and All  
Instrument Methods



## Chromatography Data Software (CDS)

Text Mix pH in S-Matrix - MD DemoLC Tutorial - Sample Workup as System/Administrator - Sample Set Method Editor


File Edit View Help

Apply Table Preferences Sample Set Method

Vial	Inj Vol (uL)	# of Injs	Label	SampleName	Level	Function	Method Set / Report Method	Label Reference	Processing	Run Time (Minutes)	Data Start (Minutes)	Next Inj. Delay (Minutes)	Column Position	Auto Additions	SampleWeight	Dilution
1						Condition Column				6.70	0.00	0.00	No Change			
2						Condition Column	Text Mix pH 001_017			0.10	0.00	0.00	No Change			
3						Equilibrate	Text Mix pH 001_017			3.00	0.00	7.95	No Change			
4	1	2.0	1	Unk-000-000	Blank - 1	Inject Samples	Text Mix pH 001_017		Normal	10.50	0.00	1.50			1.00000	1.00000
5						Condition Column	Text Mix pH 001_001			0.10	0.00	0.00	No Change			
6						Equilibrate	Text Mix pH 001_001			3.00	0.00	0.00	No Change			
7	2	2.0	1	Unk-001-001	1.a.1.a	Inject Samples	Text Mix pH 001_001		Normal	10.50	0.00	1.50			1.00000	1.00000
8						Condition Column	Text Mix pH 001_002			0.10	0.00	0.00	No Change			
9						Equilibrate	Text Mix pH 001_002			3.00	0.00	0.00	No Change			
10	2	2.0	1	Unk-001-002	2.a.1.a	Inject Samples	Text Mix pH 001_002		Normal	10.50	0.00	1.50			1.00000	1.00000
11						Condition Column				6.70	0.00	0.00	No Change			
12						Condition Column	Text Mix pH 001_003			0.10	0.00	0.00	No Change			
13						Equilibrate	Text Mix pH 001_003			3.00	0.00	0.00	No Change			
14	2	2.0	1	Unk-001-003	3.a.1.a	Inject Samples	Text Mix pH 001_003		Normal	10.50	0.00	1.50			1.00000	1.00000
15						Condition Column				6.70	0.00	0.00	No Change			
16						Condition Column	Text Mix pH 001_004			0.10	0.00	0.00	No Change			
17						Equilibrate	Text Mix pH 001_004			3.00	0.00	0.00	No Change			
18	2	2.0	1	Unk-001-004	4.a.1.a	Inject Samples	Text Mix pH 001_004		Normal	10.50	0.00	1.50			1.00000	1.00000
19						Condition Column	Text Mix pH 001_005			0.10	0.00	0.00	No Change			
20						Equilibrate	Text Mix pH 001_005			3.00	0.00	0.00	No Change			
21	2	2.0	1	Unk-001-005	5.a.1.a	Inject Samples	Text Mix pH 001_005		Normal	10.50	0.00	1.50			1.00000	1.00000
22						Condition Column				6.70	0.00	0.00	No Change			
23						Condition Column	Text Mix pH 001_006			0.10	0.00	0.00	No Change			
24						Equilibrate	Text Mix pH 001_006			3.00	0.00	0.00	No Change			

**Eliminate Transcription Errors.  
Maintain Data in Audited Environment.**

# Fusion QbD – Automated DOE Experiment Workflow




Automatically Retrieves  
All Chromatogram  
Results Data

Automated analysis,  
graphing, and reporting.

Report formats:  
RTF, DOC, HTML, PDF,  
XLSX, XML

## Chromatography Data Software (CDS)

Run	Vial	# of Inj	Label	Sample Name	Level	Function	Method Set / Report Method	Label Reference	Processing	Run Time (minutes)	Data Start (minutes)	Inlet Inj Delay (minutes)	Column Position	Auto Addition	Sample/Inj(s)	Dilution
11						Condition Column	Test file pH 81_2317			8.70	0.00	0.00	No Change			
21						Condition Column	Test file pH 81_2317			9.10	0.00	0.00	No Change			
31						Equilibrate	Test file pH 81_2317			3.00	0.00	7.00	No Change			
41	1	2.0	1	UW-001001	2.a.1.a	Injct Samples	Test file pH 81_2317		Normal	10.50	0.00	1.50	No Change		1.00000	1.00000
51						Condition Column	Test file pH 81_2317			0.10	0.00	0.00	No Change			
61						Equilibrate	Test file pH 81_2317			3.00	0.00	0.00	No Change			
71	2	2.0	1	UW-001001	1.a.1.a	Injct Samples	Test file pH 81_2317		Normal	10.50	0.00	1.50	No Change		1.00000	1.00000
81						Condition Column	Test file pH 81_2322			0.10	0.00	0.00	No Change			
91						Equilibrate	Test file pH 81_2322			3.00	0.00	0.00	No Change			
101	2	2.0	1	UW-001002	2.a.1.a	Injct Samples	Test file pH 81_2322		Normal	10.50	0.00	1.50	No Change		1.00000	1.00000
111						Condition Column	Test file pH 81_2322			8.70	0.00	0.00	No Change			
121						Condition Column	Test file pH 81_2322			0.10	0.00	0.00	No Change			
131						Equilibrate	Test file pH 81_2322			3.00	0.00	0.00	No Change			
141	2	2.0	1	UW-001001	1.a.1.a	Injct Samples	Test file pH 81_2322		Normal	10.50	0.00	1.50	No Change		1.00000	1.00000
151						Condition Column	Test file pH 81_2324			6.70	0.00	0.00	No Change			
161						Condition Column	Test file pH 81_2324			0.10	0.00	0.00	No Change			
171						Equilibrate	Test file pH 81_2324			3.00	0.00	0.00	No Change			
181	2	2.0	1	UW-001004	1.a.1.a	Injct Samples	Test file pH 81_2324		Normal	10.50	0.00	1.50	No Change		1.00000	1.00000
191						Condition Column	Test file pH 81_2325			0.10	0.00	0.00	No Change			
201						Equilibrate	Test file pH 81_2325			3.00	0.00	0.00	No Change			
211	2	2.0	1	UW-001005	1.a.1.a	Injct Samples	Test file pH 81_2325		Normal	10.50	0.00	1.50	No Change		1.00000	1.00000
221						Condition Column	Test file pH 81_2325			8.70	0.00	0.00	No Change			
231						Condition Column	Test file pH 81_2325			0.10	0.00	0.00	No Change			
241						Equilibrate	Test file pH 81_2325			3.00	0.00	0.00	No Change			



**Eliminate Transcription Errors.  
Maintain Data in Audited Environment.**